

CAB9239

Product Information

Product SKU:	CAB9239	Gene ID:	79641	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	32kDa	Conjugate:	Unconjugated
Calculated MW:	32kDa	Isotype:	IgG

Immunogen Information

Background: Involved in brain development; neurogenesis; and odontogenesis of dentin-containing tooth. Located in nuclear envelope. Implicated in Kohlschutter-Tonz syndrome.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: KTZS; RAV2; ROGDI1; ROGDI

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ROGDI (NP_078865.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.