

CAB9336

Product Information

Product SKU:	CAB9336	Gene ID:	5537	Size:	20uL, 100uL
Clone No:	ARC2709	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	34kDa	Conjugate:	Unconjugated
Calculated MW:	35kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes the catalytic subunit of protein phosphatase, a component of a signaling pathway regulating cell cycle progression. Splice variants encoding different protein isoforms exist. The pseudogene of this gene is located on chromosome X.

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: PP6; PP6C; PPP6C

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 206-305 of human PPP6C (O00743).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.