

CAB9348

Product Information

Product SKU:	CAB9348	Gene ID:	1209	Size:	20uL, 100uL
Clone No:	ARC2720	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	76kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in regulation of T cell differentiation in thymus. Located in membrane.

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: CLPTM1; transmembrane protein

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CLPTM1 (O96005).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.