

CAB9391

Product Information

Product SKU: CAB9391	Gene ID: 9918	Size: 20uL, 100uL
Clone No: ARC2757	Host Species: Rabbit	Reactivity: Human,Mouse

Additional Information

Observed MW: 157kDa	Conjugate: Unconjugated
Calculated MW: 157kDa	Isotype: IgG

Immunogen Information

Background: Enables histone binding activity. Involved in mitotic chromosome condensation. Located in condensed chromosome; cytosol; and nucleoplasm. Part of condensin complex. Colocalizes with cytoplasm and nuclear chromosome. Implicated in primary autosomal recessive microcephaly.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: CNAP1; CAP-D2; MCPH21; hCAP-D2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human CNAP1 (Q15021).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.