

CAB9434

Product Information

Product SKU:	CAB9434	Gene ID:	60491	Size:	20uL, 100uL
Clone No:	ARC2775	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background: Enables identical protein binding activity. Involved in positive regulation of transcription, DNA-templated. Located in cytoplasm.

Recommended Dilution: WB,1:100 - 1:500

Synonyms: CALS-7; MDS015; ALS2CR1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 250-350 of human ALS2CR1 (Q9GZT8).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.