

## CAB9546

---

### Product Information

<b>Product SKU:</b>	CAB9546	<b>Gene ID:</b>	1044	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC1627	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	28kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	28kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded DNA-binding protein regulates intestine-specific gene expression and enterocyte differentiation. It has been shown to induce expression of the intestinal alkaline phosphatase gene, and inhibit beta-catenin/T-cell factor transcriptional activity.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** -

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cdx1 (P47902).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.