

## CAB9580

### Product Information

<b>Product SKU:</b>	CAB9580	<b>Gene ID:</b>	6434	<b>Size:</b>	20uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

### Additional Information

<b>Observed MW:</b>	38kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	34kDa	<b>Isotype:</b>	IgG

### Immunogen Information

**Background:** This gene encodes a nuclear protein which functions as sequence-specific serine/arginine splicing factor which plays a role in mRNA processing, splicing patterns, and gene expression. Alternative splicing results in multiple transcript variants.

**Recommended Dilution:** WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

**Synonyms:** SFRS10; SRFS10; TRAN2B; PPP1R156; TRA2-BETA; Htra2-beta; TRA2B

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide of human TRA2B

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.