

CAB9701

Product Information

Product SKU:	CAB9701	Gene ID:	10928	Size:	20uL, 100uL
Clone No:	ARC1707	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	-	Conjugate:	Unconjugated
Calculated MW:	76kDa	Isotype:	IgG

Immunogen Information

Background:	RALBP1 plays a role in receptor-mediated endocytosis and is a downstream effector of the small GTP-binding protein RAL (see RALA; MIM 179550). Small G proteins, such as RAL, have GDP-bound inactive and GTP-bound active forms, which shift from the inactive to the active state through the action of RALGDS (MIM 601619), which in turn is activated by RAS (see HRAS; MIM 190020) (summary by Feig, 2003 [PubMed 12888294]).
Recommended Dilution:	IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200
Synonyms:	RIP1; RLIP1; RLIP76; RALBP1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human RALBP1 (Q15311).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.