

CAB9772

Product Information

Product SKU: CAB9772	Gene ID: 997	Size: 20uL, 100uL
Clone No: ARC1742	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 32kDa	Conjugate: Unconjugated
Calculated MW: 27kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is a member of the ubiquitin-conjugating enzyme family. Ubiquitin-conjugating enzyme catalyzes the covalent attachment of ubiquitin to other proteins. This protein is a part of the large multiprotein complex, which is required for ubiquitin-mediated degradation of cell cycle G1 regulators, and for the initiation of DNA replication.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: UBC3; UBCH3; UBE2R1; E2-CDC34; CDC34

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-236 of human CDC34 (P49427).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.