

## CABC022

---

### Product Information

<b>Product SKU:</b>	CABC022	<b>Gene ID:</b>	-	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Various Plants

---

### Additional Information

<b>Observed MW:</b>	42kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	45kDa	<b>Isotype:</b>	IgG1

---

### Immunogen Information

**Background:** Actins are highly conserved proteins that are involved in cell motility, structure, and integrity. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** -

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to Actin (plant specific).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.