

CABC042

Product Information

Product SKU:	CABC042	Gene ID:	-	Size:	20uL, 100uL
Clone No:	ARC5105-03	Host Species:	Rabbit	Reactivity:	-

Additional Information

Observed MW:	-	Conjugate:	Unconjugated
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a transcriptional regulator and tumor suppressor, serving as an activator of genes involved in both innate and acquired immune responses.

Recommended Dilution: IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200 FC,1:50 - 1:200 IP,0.5µg-4µg antibody for 200µg-400µg extracts of whole cellsChIP,3µg antibody for 10µg-15µg of Chromatin

Synonyms: -

Purification Method: Protein A/G purification

Immunogen: None

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.