

CABC053

Product Information

Product SKU:	CABC053	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Widely Range, Oryza sativa

Additional Information

Observed MW:	45kDa	Conjugate:	Unconjugated
Calculated MW:	45kDa	Isotype:	IgG1

Immunogen Information

Background: Actins are highly conserved proteins that are involved in cell motility, structure, and integrity. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins.

Recommended Dilution: WB, 1:1000 - 1:5000

Synonyms: ACT12; Actin(Plant Specific)

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of Actin(Plant Specific)(NP_190236.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.