

Phospho-MAPK14-T180/Y182+MAPK11-T180/Y182+MAPK13-T180/Y182+MAPK12-T183/Y1845 Rabbit Polyclonal Antibody CABP1256

Product Information

Product SKU:	CABP1256	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	43kDa	Conjugate:	Unconjugated
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background: Mapk14 (phosphoT180/Y182) monoclonal antibody is provided, clone number M139 is a monoclonal antibody prepared from mice, its clone number is M139. Mapk14 is also known as CSBP, CSPB1, Csbp1, Csbp2, Exip, Hog, MGC105413, Mxi2, Prkm14, Prkm15, RK, Sapk2A, p38, p38Hog, p38alpha; The full name of its gene is Mitogen Activated Protein Kinase 14.

Recommended Dilution: WB,1:100 - 1:500

Synonyms: -

Purification Method: Affinity purification

Immunogen: A synthetic phosphorylated peptide around T180 & Y182 of human MAPK14MAPK11 (NP_002742.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.