

---

## Product Information

<b>Product SKU:</b> CABP1284	<b>Gene ID:</b> 7430/5962/4478	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human

---

## Additional Information

<b>Observed MW:</b> 75kDa/80kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 69kDa/68kDa/67kDa	<b>Isotype:</b> IgG

---

## Immunogen Information

**Background:** The cytoplasmic peripheral membrane protein encoded by this gene functions as a protein-tyrosine kinase substrate in microvilli. As a member of the ERM protein family, this protein serves as an intermediate between the plasma membrane and the actin cytoskeleton. This protein plays a key role in cell surface structure adhesion, migration and organization, and it has been implicated in various human cancers. A pseudogene located on chromosome 3 has been identified for this gene. Alternatively spliced variants have also been described for this gene. Radixin is a cytoskeletal protein that may be important in linking actin to the plasma membrane. It is highly similar in sequence to both ezrin and moesin. The radixin gene has been localized by fluorescence in situ hybridization to 11q23. A truncated version representing a pseudogene (RDXP2) was assigned to Xp21.3. Another pseudogene that seemed to lack introns (RDXP1) was mapped to 11p by Southern and PCR analyses. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. Moesin (for membrane-organizing extension spike protein) is a member of the ERM family which includes ezrin and radixin. ERM proteins appear to function as cross-linkers between plasma membranes and actin-based cytoskeletons. Moesin is localized to filopodia and other membranous protrusions that are important for cell-cell recognition and signaling and for cell movement.

<b>Recommended Dilution:</b>	WB, 1:500 - 1:1000
<b>Synonyms:</b>	CVL; CVIL; VIL2; HEL-S-105; DFNB24; HEL70; IMD50; Phospho-Ezrin-Thr567/Radixin-Thr564/Moesin-Thr558
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	A synthetic phosphorylated peptide around T567 of human Ezrin/Radixin/MoesinEZR (NP_003370.2).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.