

## CABP1318

---

### Product Information

<b>Product SKU:</b>	CABP1318	<b>Gene ID:</b>	5970	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	65kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	60kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** NF-kappa-B is a ubiquitous transcription factor involved in several biological processes. It is held in the cytoplasm in an inactive state by specific inhibitors. Upon degradation of the inhibitor, NF-kappa-B moves to the nucleus and activates transcription of specific genes. NF-kappa-B is composed of NFKB1 or NFKB2 bound to either REL, RELA, or RELB. The most abundant form of NF-kappa-B is NFKB1 complexed with the product of this gene, RELA. Four transcript variants encoding different isoforms have been found for this gene.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** p65; CMCU; NFKB3; AIF3BL3; Phospho-NF-kB p65/RelA-T505

**Purification Method:** Affinity purification

**Immunogen:** A synthetic phosphorylated peptide around T505 of human Phospho-NF-kB p65/RelA-T505 (NP\_068810.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.