

CABP1324

Product Information

Product SKU:	CABP1324	Gene ID:	4137	Size:	20uL, 100uL
Clone No:	ARC54951	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	50-80kDa	Conjugate:	Unconjugated
Calculated MW:	79kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes the microtubule-associated protein tau (MAPT) whose transcript undergoes complex, regulated alternative splicing, giving rise to several mRNA species. MAPT transcripts are differentially expressed in the nervous system, depending on stage of neuronal maturation and neuron type. MAPT gene mutations have been associated with several neurodegenerative disorders such as Alzheimer's disease, Pick's disease, frontotemporal dementia, cortico-basal degeneration and progressive supranuclear palsy.

Recommended Dilution: WB,1:2000 - 1:20000

Synonyms: TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T181

Purification Method: Affinity purification

Immunogen: A synthetic phosphorylated peptide around T181 of human Phospho-Tau-T181 (NP_001190181.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.