

CABP1355

Product Information

Product SKU:	CABP1355	Gene ID:	4790	Size:	20uL, 100uL
Clone No:	ARC55043	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	120kDa	Conjugate:	Unconjugated
Calculated MW:	105kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a 105 kD protein which can undergo cotranslational processing by the 26S proteasome to produce a 50 kD protein. The 105 kD protein is a Rel protein-specific transcription inhibitor and the 50 kD protein is a DNA binding subunit of the NF-kappa-B (NFKB) protein complex. NFKB is a transcription regulator that is activated by various intra- and extra-cellular stimuli such as cytokines, oxidant-free radicals, ultraviolet irradiation, and bacterial or viral products. Activated NFKB translocates into the nucleus and stimulates the expression of genes involved in a wide variety of biological functions. Inappropriate activation of NFKB has been associated with a number of inflammatory diseases while persistent inhibition of NFKB leads to inappropriate immune cell development or delayed cell growth. NFKB is a critical regulator of the immediate-early response to viral infection. Alternative splicing results in multiple transcript variants encoding different isoforms, at least one of which is proteolytically processed.

Recommended Dilution: WB,1:1000 - 1:5000

Synonyms: KBF1; EBP-1; NF-kB; CVID12; NF-kB1; NFKB-p50; NFkappaB; NF-kappaB; NFKB-p105; NF-kappa-B1; NF-kappabeta; Phospho-NF- κ B p105-S932

Purification Method: Affinity purification

Immunogen: A synthetic phosphorylated peptide around S932 of human Phospho-NF- κ B p105-S932 (NP_001158884.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.