

CABP1411

Product Information

Product SKU:	CABP1411	Gene ID:	4090	Size:	20uL, 100uL
Clone No:	ARC59204	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	58kDa	Conjugate:	Unconjugated
Calculated MW:	52kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is involved in the transforming growth factor beta signaling pathway that results in an inhibition of the proliferation of hematopoietic progenitor cells. The encoded protein is activated by bone morphogenetic proteins type 1 receptor kinase, and may be involved in cancer. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB,1:1000 - 1:5000

Synonyms: DWFC; JV5-1; MADH5; Phospho-Smad5-S465

Purification Method: Affinity purification

Immunogen: A synthetic phosphorylated peptide around S465 of human Smad5.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.