

CABP1432

Product Information

Product SKU:	CABP1432	Gene ID:	3991	Size:	20uL, 100uL
Clone No:	ARC61254	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	84kDa	Conjugate:	Unconjugated
Calculated MW:	117kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene has a long and a short form, generated by use of alternative translational start codons. The long form is expressed in steroidogenic tissues such as testis, where it converts cholesteryl esters to free cholesterol for steroid hormone production. The short form is expressed in adipose tissue, among others, where it hydrolyzes stored triglycerides to free fatty acids.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: HSL; LHS; REH; AOMS4; FPLD6; Phospho-HSL-S660

Purification Method: Affinity purification

Immunogen: A synthesized peptide derived from human HSL (Phospho-Ser660)

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH 7.3.