

CABP1510

Product Information

Product SKU:	CABP1510	Gene ID:	7408	Size:	20uL, 100uL
Clone No:	ARC64203	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	50kDa	Conjugate:	Unconjugated
Calculated MW:	40kDa	Isotype:	IgG

Immunogen Information

Background:	Vasodilator-stimulated phosphoprotein (VASP) is a member of the Ena-VASP protein family. Ena-VASP family members contain an EHV1 N-terminal domain that binds proteins containing E/DFPPPPXD/E motifs and targets Ena-VASP proteins to focal adhesions. In the mid-region of the protein, family members have a proline-rich domain that binds SH3 and WW domain-containing proteins. Their C-terminal EVH2 domain mediates tetramerization and binds both G and F actin. VASP is associated with filamentous actin formation and likely plays a widespread role in cell adhesion and motility. VASP may also be involved in the intracellular signaling pathways that regulate integrin-extracellular matrix interactions. VASP is regulated by the cyclic nucleotide-dependent kinases PKA and PKG.
Recommended Dilution:	WB,1:2000 - 1:5000
Synonyms:	VASP; Phospho-VASP-S157
Purification Method:	Affinity purification
Immunogen:	A synthetic phosphorylated peptide around S157 of human VASP(NP_003361.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.