

CABP1513

Product Information

Product SKU:	CABP1513	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Other (Wide Range Predicted)

Additional Information

Observed MW:	10-250kDa/>10kDa	Conjugate:	Unconjugated
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background:

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence containing phosphorylated S & T.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.