

Dengue type 2 Antibody



CPAB0381

Product Information

Size:

500µg

Applications:**Reactivity:**

Viral

Source:

Mouse

Isotype:

IgG1

Purification Method:

Ion exchange column.

Protein Background:

Caused by one of four closely related virusserotypesof the genus Flavivirus, family Flaviviridae, each serotype is sufficiently different that there is no cross-protection and epidemics caused by multiple serotypes (hyperendemicity) can occur. In cell culture experiments and mice Morpholinoantisense oligos have shown specific activity against Dengue virus.

Synonyms:**Immunogen:**

Dengue-infected cell lysate

Storage:

In lyophilized form, for long periods, store at 4ºC in a dry environment. After reconstitution, if not intended for use within a month, aliquot and store at -20ºC.