

Influenza-A H5N1 Antibody



CPAB0390

Product Information

Size:

10µg

Applications:

Western Blot, ELISA

Reactivity:

Viral

Source:**Isotype:****Purification Method:**

Chromatography on protein G
Sepharose.

Protein Background:

Influenza A virus subtype H5N1 is a subtype of the Influenza A virus. Subtype H5N1 is commonly known as avian influenza or bird flu. H5N1 may cause more than one influenza pandemic as it is expected to continue mutating in birds. The dominant strain of HPAI A (H5N1) evolved creating the Z genotype. It has also been called Asian lineage HPAI/A/H5N1 which is divided into 2 antigenic clades. "Clade 1 includes human and bird isolates from Vietnam, Thailand, and Cambodia and bird isolates from Laos and Malaysia. Clade 2 viruses were first identified in bird isolates from China, Indonesia, Japan, and South Korea before spreading westward to the Middle East, Europe, and Africa.

Synonyms:**Immunogen:**

Influenza A hemagglutinin H5N1.

Storage:

Should be stored at 4C. Do not freeze.