

Dengue NS1 Antibody



CPAB0396

Product Information

Size:

100µg

Applications:**Reactivity:**

Viral

Source:**Isotype:****Purification Method:**

Purified monoclonal IgG by protein A chromatography.

Protein Background:

Paired monoclonal antibodies developed for dengue NS1 antigen test by lateral flow immunoassay device. First antibody is for conjugation and the second antibody is for membrane coating. Please note that when ordering for example: 100µg antibody we ship 50µg from each of the antibody solutions (100µg in total).

Synonyms:**Immunogen:****Storage:**