

Dengue NS1 Paired Antibody



CPAB0397

Product Information

Size:

100µg

Applications:

Reactivity:

Viral

Source:

Isotype:

Protein Background:

Dengue NS1 is an early marker for dengue infection from the first day of fever and lasting about one week. Dengue IgM is usually detected early on the first day to 5 days of infection. Dengue NS1 could be detected from both primary and secondary infection for all four subtypes.

Synonyms:

Immunogen:

Storage:

Purification Method:

Purified monoclonal IgG by protein A chromatography.