

# LY6C/G Antibody, Biotin



CPAB0454

---

## Product Information

**Size:**

500µg

**Applications:**

Immunohistochemistry, Flow Cytometry

**Reactivity:**

Other bodies

**Source:**

Rat

**Isotype:**

IgG2b

**Purification Method:**

Protein A column.

**Protein Background:**

**Synonyms:**

**Immunogen:**

BM-enriched myeloid cells.

**Storage:**

In lyophilized form, for long periods, store at 4°C in a dry environment. After reconstitution, if not intended for use within a month, aliquot and store at -20°C.