

Chlamydia antibody



CPAB0464

Product Information

Size:

500µg

Applications:

Reactivity:

Other bodies

Source:

Mouse

Isotype:

IgG2b

Purification Method:

Ion exchange column.

Protein Background:

Chlamydia is a common term for infection with any bacterium belonging to the phylum Chlamydiae. This term derives from the name of the bacterial genus Chlamydia in the family Chlamydiaceae, order Chlamydiales, class and phylum Chlamydiae. There are two genera in Chlamydiaceae: Chlamydia and Chlamydophila. The genus Chlamydia includes three species: C. trachomatis, C. muridarum, and C. suis.

Synonyms:

Immunogen:

Purified Chlamydia LPS.

Storage:

In lyophilized form, for long periods, store at 4°C in a dry environment. After reconstitution, if not intended for use within a month, aliquot and store at -20°C.