

# Dengue 3 NS1 Antibody



CPAB0704

---

## Product Information

**Size:**

100µg

**Applications:**

**Reactivity:**

Viral

**Source:**

**Isotype:**

**Purification Method:**

Purified by protein A chromatography.

**Protein Background:**

The polyclonal antibody to dengue serotype 3 NS1 was collected from the rabbit immunized with full length recombinant dengue serotype 3 NS1 antigen, reactive to all dengue serotype NS1.

**Synonyms:**

**Immunogen:**

**Storage:**