

Anti-UCHL1 Biotinylated Recombinant Rabbit Monoclonal Antibody

Product Information

Product SKU:

HDAB0216

Size:

100 µg

Isotype:

Rabbit IgG

Host:

Rabbit

Clonality:

Monoclonal

Clone:

DM214

Category:

Recombinant Antibody

Reactivity:

Human

Synonyms:

HEL-117; HEL-S-53; NDGOA; PARK5; PGP 9.5; PGP9.5; PGP95; SPG79; Uch-L1

Formulation:

Powder

Purification Method:

Purified from cell culture supernatant by affinity chromatography

Buffer:

1XPBS

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

ELISA, Flow Cyt

Recommended Dilutions:

ELISA 1:5000-10000; Flow Cyt 1:100

Background:

The protein encoded by this gene belongs to the peptidase C12 family. This enzyme is a thiol protease that hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin. This gene is specifically expressed in the neurons and in cells of the diffuse neuroendocrine system. Mutations in this gene may be associated with Parkinson disease. [provided by RefSeq, Sep 2009]

Research Use Only:

Research Use Only