

Anti-CD96 Chimeric Recombinant Rabbit Monoclonal Antibody

Product Information

Product SKU:

HDAB0243

Size:

100 µg

Isotype:

Rabbit/Human Fc chimeric IgG1

Host:

Rabbit

Clonality:

Monoclonal

Clone:

DMC285

Category:

Recombinant Antibody

Reactivity:

Human

Synonyms:

TACTILE

Formulation:

Powder

Purification Method:

Purified from cell culture supernatant by affinity chromatography

Buffer:

1XPBS

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

Flow Cyt

Recommended Dilutions:

Flow Cyt 1:100

Background:

The protein encoded by this gene belongs to the immunoglobulin superfamily. It is a type I membrane protein. The protein may play a role in the adhesive interactions of activated T and NK cells during the late phase of the immune response. It may also function in antigen presentation. Alternative splicing generates multiple transcript variants encoding distinct isoforms.

Research Use Only:

Research Use Only