

# Anti-SLAMF5 Chimeric Recombinant Rabbit Monoclonal Antibody

## Product Information

**Product SKU:**

HDAB0244

**Size:**

100 µg

**Isotype:**

Rabbit/Human Fc chimeric IgG1

**Host:**

Rabbit

**Clonality:**

Monoclonal

**Clone:**

DMC286

**Category:**

Recombinant Antibody

**Reactivity:**

Human

**Synonyms:**

CD84,SLAMF5,LY9B,SLAMF5

**Formulation:**

Powder

**Purification Method:**

Purified from cell culture supernatant by affinity chromatography

**Buffer:**

1XPBS

**Storage and Stability:**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

**Applications**

Flow Cyt

**Recommended Dilutions:**

Flow Cyt 1:100

**Background:**

This gene encodes a membrane glycoprotein that is a member of the signaling lymphocyte activation molecule (SLAM) family. This family forms a subset of the larger CD2 cell-surface receptor Ig superfamily. The encoded protein is a homophilic adhesion molecule that is expressed in numerous immune cells types and is involved in regulating receptor-mediated signaling in those cells. Alternate splicing results in multiple transcript variants.

**Research Use Only:**

Research Use Only