

# Anti-SLAMF1 Chimeric Recombinant Rabbit Monoclonal Antibody

## Product Information

**Product SKU:**

HDAB0283

**Size:**

100 µg

**Isotype:**

Rabbit/Human Fc chimeric IgG1

**Host:**

Rabbit

**Clonality:**

Monoclonal

**Clone:**

DMC445

**Category:**

Recombinant Antibody

**Reactivity:**

Human

**Synonyms:**

SLAMF1, CD150, CDw150, SLAM, IPO-3

**Formulation:**

Powder

**Purification Method:**

Purified from cell culture supernatant by affinity chromatography

**Buffer:**

1XPBS

**Storage and Stability:**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

**Applications**

Flow Cyt

**Recommended Dilutions:**

Flow Cyt 1:100

**Background:**

SLAMF1 (Signaling Lymphocytic Activation Molecule Family Member 1) is a Protein Coding gene. Diseases associated with SLAMF1 include Measles and Subacute Sclerosing Panencephalitis. Among its related pathways are Dendritic Cells Developmental Lineage Pathway and T Cell Co-Signaling Pathway: Ligand-Receptor Interactions. Gene Ontology (GO) annotations related to this gene include antigen binding. An important paralog of this gene is LY9.

**Research Use Only:**

Research Use Only