

Anti-CRTAM Chimeric Recombinant Rabbit Monoclonal Antibody

Product Information

Product SKU:

HDAB0324

Size:

100 µg

Isotype:

Rabbit/Human Fc chimeric IgG1

Host:

Rabbit

Clonality:

Monoclonal

Clone:

DMC673

Category:

Recombinant Antibody

Reactivity:

Human

Synonyms:

CD355

Formulation:

Powder

Purification Method:

Purified from cell culture supernatant by affinity chromatography

Buffer:

1XPBS

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

Flow Cyt

Recommended Dilutions:

Flow Cyt 1:100

Background:

The CRTAM gene is upregulated in CD4 (see MIM 186940)-positive and CD8 (see CD8A; MIM 186910)-positive T cells and encodes a type I transmembrane protein with V and C1-like Ig domains (Yeh et al., 2008 [PubMed 18329370]). [supplied by OMIM, Feb 2009]

Research Use Only:

Research Use Only

AssayGenie

