

Tezepelumab Biosimilar (Anti-TSLP) Antibody

Product Information

Product SKU:

HDBS0111

Size:

100 µg

Isotype:

IgG2

Host:

Homo sapiens

Category:

Biosimilar Antibody

Uniprot:

Q969D9

Reactivity:

Human

Synonyms:

Thymic stromal lymphopoietin

Common Name:

AMG 157, AMG-157, MEDI-9929, MEDI9929, anti-TSLP-Receptor (TSLP-R; CrI2), tezepelumab-ekko

Target:

TSLP

Formulation:

Powder

Buffer:

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

Guarantee:

12 months from date of dispatch

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

ELISA, Flow Cyt

Recommended Dilutions:

ELISA 1:5000-10000; Flow Cyt 1:100

Research Use Only:

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.