

## Lirilumab Biosimilar (Anti-KIR2DL2) Antibody

### Product Information

**Product SKU:**

HDBS0121

**Size:**

100 µg

**Isotype:**

IgG4

**Host:**

Homo sapiens

**Category:**

Biosimilar Antibody

**Uniprot:**

P43627

**Reactivity:**

Human

**Synonyms:**

CD158b, CD158B1, NKAT-6, NKAT6, p58.2

**Common Name:**

BMS-986015, IPH2102

**Target:**

KIR2DL2

**Formulation:**

Powder

**Buffer:**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.



**Guarantee:**

12 months from date of dispatch

**Storage and Stability:**

Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

**Applications**

Flow Cyt

**Recommended Dilutions:**

Flow Cyt 1:100

**Research Use Only:**

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

