

Vixarelimab Biosimilar (Anti-OSMR) Antibody

Product Information

Product SKU:

HDBS0126

Size:

100 µg

Isotype:

IgG4-lambda2

Host:

Homo sapiens

Category:

Biosimilar Antibody

Uniprot:

Q99650

Reactivity:

Human

Synonyms:

IL-31 receptor subunit beta, IL-31R subunit beta, IL-31R-beta, IL-31RB

Common Name:

KPL-716

Target:

OSMR

Formulation:

Powder

Buffer:

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.



Guarantee:

12 months from date of dispatch

Storage and Stability:

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Applications

Flow Cyt

Recommended Dilutions:

Flow Cyt 1:100

Research Use Only:

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

