

## IVMB0457

---

### Product Information

<b>Product SKU:</b> IVMB0457	<b>Clone:</b> Hu137	<b>Target:</b> IL-6
<b>Size:</b> 100 µg		<b>Isotype:</b> Human IgG1κ

---

### Additional Information

<b>Reactivity:</b> Human	<b>Host Species:</b> Human
<b>Antibody Type:</b> Biosimilar Recombinant Human Monoclonal Antibody	<b>Expression Host:</b> HEK-293 Cells

---

### Immunogen Information

**Background:** IL-6 and its signaling pathway play a part in immune response regulation, inflammation, and hematopoiesis.<sup>2</sup> Sarilumab is a research-grade recombinant human monoclonal IL-6 receptor antagonist. It specifically binds to both the transmembrane and soluble forms of the IL-6 receptor, thus inhibiting IL-6-mediated cis and trans-signaling in a dose-dependent manner.<sup>1</sup> Therapeutic Sarilumab, also known by the trade name Kevzara, is currently used to treat Rheumatoid Arthritis<sup>1</sup>, however, as of March 2020, The Feinstein Institute of Northwell Health publicized a study on "a human antibody that may prevent the activity" of IL-6 for the treatment of COVID-19.<sup>3</sup> Anti-Human IL-6 (Sarilumab) utilizes the same variable regions from the therapeutic antibody Sarilumab making it ideal for research projects.

**Product Concentration:** 0.5 mg/ml

**Applications:** ELISA

**Synonyms:** Interleukin-6, CDF; HGF; HSF; BSF2; BSF-2; IFNB2; IFN-beta-2

**Antigen Distribution:** IL-6R is ubiquitously expressed.

**Additional Applications:** FA

**Formulation:** This Biotinylated antibody is formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.4, 1% BSA and 0.09% sodium azide as a preservative.

**Contact Details | Dublin, Ireland**

**Email:** [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | **Web:** [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.

- Specificity:** This non-therapeutic antibody uses the same variable region sequence as the therapeutic antibody Sarilumab. Sarilumab binds to the transmembrane and soluble forms of the IL-6 receptor. This product is for research use only.
- Pathogen Testing:** -
- Storage & Handling:** This biotinylated antibody is stable when stored at 2-8°C. Do not freeze.