

# CD1A Monoclonal Antibody



MAC00042

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Mouse

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:200

**Protein Background:**

**Gene ID:**

CD1A

**Uniprot**

P06126

**Synonyms:**

T-cell surface glycoprotein CD1a

**Immunogen:**

Synthetic Peptide

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---

N/A

N/A