

Tri-methyl-Histone H3(K9) Monoclonal Antibody



MAC0071

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:**Gene ID:**

Tri-methyl-Histone H3(K9)

Uniprot

P68431/Q71DI3/P84243

Synonyms:

HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G;


Immunogen:

Synthetic Peptide

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) Hela, 2) Rat Testis tissue, 3) Raw264.7, diluted at 1:2000.