

# Di-methyl-Histone H3(K36) Monoclonal Antibody



MACO0087

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Mouse

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:5000

**Protein Background:****Gene ID:**

Di-methyl-Histone H3(K36)

**Uniprot**

P68431/Q71DI3/P84243

**Synonyms:**

HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G;

**Immunogen:**


Synthetic Peptide

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of 1) HeLa, 2) Raw264.7, 3) Rat Brain Tissue, diluted at 1:3,000.