

# TUBB3 Monoclonal Antibody



MAC00122


---

## Product Information

<b>Size:</b> 50ug	<b>Protein Background:</b>
<b>Reactivity:</b> Human, Mouse, Rat	<b>Gene ID:</b> TUBB3
<b>Source:</b> Mouse	<b>Uniprot</b> Q13509
<b>Isotype:</b> IgG	<b>Synonyms:</b> TUBB3; TUBB4; Tubulin beta-3 chain; Tubulin beta-4 chain; Tubulin beta-III
<b>Applications:</b> ELISA, WB	<b>Immunogen:</b> Synthetic Peptide
<b>Recommended dilutions:</b> WB:1:500-1:5000	<b>Storage:</b> PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of 1) HeLa, 2) Mouse Brain, 3) Rat Brain, diluted at 1:5000.