

EBFP Monoclonal Antibody



MAC00128

Product Information

Size:

50ug

Protein Background:

Reactivity:

Gene ID:

EBFP

Source:

Mouse

Uniprot

Isotype:

IgG

Synonyms:

Applications:

ELISA, WB

Immunogen:

Recombinant Protein

Recommended dilutions:

WB:1:500-1:5000

Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of recombinant EBFP protein, diluted at 1:5000.