

CHD3 Monoclonal Antibody



MACO0203

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Mouse

Isotype:

IgG1

Applications:

ELISA, WB, ICC

Recommended dilutions:

WB:1:500-1:5000, ICC:1:50-1:200

Protein Background:

Component of the histone deacetylase NuRD complex which participates in the remodeling of chromatin by deacetylating histones. Required for anchoring centrosomal pericentrin in both interphase and mitosis, for spindle organization and centrosome integrity.

Gene ID:

CHD3

Uniprot

Q12873

Synonyms:

ZFH; Mi-2a; Mi2-ALPHA

Immunogen:

Purified recombinant human CHD3 protein fragments expressed in E.coli.

Storage:

Ascites

Product Images

N/A

N/A