

ITGB3 Monoclonal Antibody



MAC00246

Product Information

Size:

100ul

Reactivity:

Human, Dog, Rabbit, Rat, Mouse

Source:

Mouse

Isotype:

IgG1

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:

Integrin alpha-V/beta-3 is a receptor for cytotactin, fibronectin, laminin, matrix metalloproteinase-2, osteopontin, osteomodulin, prothrombin, thrombospondin, vitronectin and von Willebrand factor. Integrin alpha-IIb/beta-3 is a receptor for fibronectin.

Gene ID:

ITGB3

Uniprot

P05106

Synonyms:

GT, CD61, GP3A, BDPLT2, GPIIIa

Immunogen:

Purified recombinant human ITGB3 protein fragments expressed in E.coli.

Storage:

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

Product Images

N/A

N/A