

# CD4/CD25 Monoclonal Antibody



MAC00295

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Mouse

**Isotype:**

IgG2b/IgG2a

**Applications:**

ELISA, IF, FC

**Recommended dilutions:**

IF:1:100-1:200

**Protein Background:**

**Gene ID:**

CD4/CD25

**Uniprot**

P01730

**Synonyms:**

**Immunogen:**

Recombinant Protein

**Storage:**

Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA

## Product Images

---

N/A

N/A