

NCAM1 Monoclonal Antibody



MACO0305

Product Information

Size:

100ul

Reactivity:

Human

Source:

Mouse

Isotype:

IgG1

Applications:

ELISA, IF, FC

Recommended dilutions:

IF:1:100-1:200

Protein Background:

COC56 reacts with CD56, a 175-220 kDa Neural Cell Adhesion Molecule (NCAM), expressed on 10-25% of peripheral blood lymphocytes, including all CD16+ NK cells and approximately 5% of CD3+ lymphocytes, referred to as NKT cells. It also is present at brain and neuromuscular junctions, certain LGL leukemias, small cell lung carcinomas, neuronally derived tumors, myeloma and myeloid leukemias. CD56 (NCAM) is involved in neuronal homotypic cell adhesion which is implicated in neural development, and in cell differentiation during embryogenesis.

Gene ID:

NCAM1

Uniprot

P13591

Synonyms:**Immunogen:**

Recombinant Protein

Storage:

Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA

Product Images

N/A

N/A