

IgG2a Monoclonal Antibody



MAC00321

Product Information

Size:

100ul

Reactivity:

Human

Source:

Mouse

Isotype:

IgG2a

Applications:

ELISA, IF, FC

Recommended dilutions:

IF:1:100-1:200

Protein Background:

This antibody can be used as a mouse IgG2a isotype control in flow cytometry and other applications.

Gene ID:

IgG2a

Uniprot

Synonyms:

Immunogen:

Recombinant Protein

Storage:

Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA

Product Images

N/A

N/A