

# PPP1R18 Monoclonal Antibody



MACO0333

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Mouse

**Isotype:**

IgG1

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:5000

**Protein Background:**

May target protein phosphatase 1 to F-actin cytoskeleton.

**Gene ID:**

PPP1R18

**Uniprot**

Q6NYC8

**Synonyms:**

HKMT1098; KIAA1949; Protein phosphatase 1 F actin cytoskeleton targeting subunit; PPP1R18;

**Immunogen:**

Purified recombinant human PPP1R18 protein fragments expressed in E.coli.

**Storage:**

Ascites

## Product Images

---

N/A

N/A