

# TBL1Y Monoclonal Antibody



MACO0363

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Pig, Rat, Dog

**Source:**

Mouse

**Isotype:**

IgG2b

**Applications:**

ELISA, WB, ICC

**Recommended dilutions:**

WB:1:500-1:5000, ICC:1:50-1:200

**Protein Background:**

**Gene ID:**

TBL1Y

**Uniprot**

Q9BQ87

**Synonyms:**

F-box-like/WD repeat-containing protein TBL1Y transducin beta-like 1 transducin beta-like protein 1Y transducin-beta-like protein 1, Y-linked

**Immunogen:**

Recombinant full length protein (522aa) of human TBL1Y.

**Storage:**

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

## Product Images

---

N/A

N/A