

Nano-Tag9 Monoclonal Antibody



MAC00375

Product Information

Size:

50ug

Reactivity:

Source:

Mouse

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:

Gene ID:

Nano-Tag9

Uniprot

Synonyms:

Immunogen:

Synthetic Peptide

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of Recombinant Nano-Tag9 Protein with Nano-Tag9 Mouse mAb diluted at 1)1:10,000 2)1:5,000